

















# Ablation Zone Chart

Disposable Microwave Therapeutic Antennas






# Ex Vivo Porcine Liver Ablation Data 14G Ceramic Antenna on Continuous Mode (for lesion<40mm)

	30W(mm)	40W(mm)	60W(mm)	80W(mm)
1min				
3min				
5min			<p><b>Note:</b> Output power 50W~80W is routinely used in operation, and the continuous working time of a single ablation circle shall not be more than 15min; The maximum power limit is 80W. It is <b>FORBIDDEN to overpower the antenna.</b></p> <p>Single antenna liver ablation is generally suitable for lesions <b>within 4cm</b>; If the lesion is larger than 4cm, double antenna ablation is recommended.</p>	
8min				
10min				
12min				





Liver








Ex Vivo Porcine Liver Ablation Data  
14G Ceramic Antenna on Pulse Mode (for lesion < 40mm)

	30W(mm)	40W(mm)	60W(mm)	80W(mm)
3min	<p>15mm 25mm</p>	<p>19mm 26mm</p>	<p>21mm 35mm</p>	<p>23mm 38mm</p>
5min	<p>20mm 27mm</p>	<p>26mm 30mm</p>	<p>27mm 37mm</p>	<p>28mm 43mm</p>
8min	<p>26mm 30mm</p>	<p>29mm 35mm</p>	<p>32mm 40mm</p>	<p>39mm 48mm</p>
10min	<p>29mm 34mm</p>	<p>30mm 40mm</p>	<p>35mm 46mm</p>	
12min	<p>30mm 35mm</p>	<p>32mm 43mm</p>		





Liver



## Note:

Output power 50W~80W is routinely used in operation, and the continuous working time of a single ablation circle shall not be more than 15min;The maximum power limit is 80W.

















It is **FORBIDDEN** to overpower the antenna.

Single antenna liver ablation is generally suitable for lesions **within 4cm**; If the lesion is larger than 4cm, double antenna ablation is recommended.






Ex Vivo Porcine Liver Ablation Data  
16G Ceramic Antenna on Continuous Mode (for lesion < 30mm)

	20W(mm)	30W(mm)	40W(mm)	50W(mm)
1min				
3min				
5min				
8min				
10min				





Liver



## Note:

Output power 30W~40W is routinely used in operation, and the continuous working time of a single ablation circle shall not be more than 15min;The maximum power limit is 50W.

**It is FORBIDDEN to overpower the antenna.**

Single antenna liver ablation is generally suitable for lesions **within 3cm**; If the lesion is larger than 3cm, double antenna ablation is recommended.





# Ex Vivo Porcine Liver Ablation Data 16G Ceramic Antenna on Pulse Mode



Liver

	20W(mm)	30W(mm)	40W(mm)	50W(mm)
1min	<p>9mm 15mm</p>	<p>10mm 20mm</p>	<p>12.5mm 24mm</p>	<p>14mm 25mm</p>
3min	<p>12mm 20mm</p>	<p>16mm 24mm</p>	<p>20mm 30mm</p>	<p>21mm 31mm</p>
5min	<p>16mm 22mm</p>	<p>18mm 27mm</p>	<p>22mm 34mm</p>	<p>23mm 35mm</p>
8min	<p>21mm 25mm</p>	<p>22mm 31mm</p>	<p>31mm 37mm</p>	<p>33mm 38mm</p>
10min	<p>22mm 26mm</p>	<p>24mm 35mm</p>	<p>32mm 43mm</p>	<p>35mm 45mm</p>

## Note:

Output power 30W~40W is routinely used in operation, and the continuous working time of a single ablation circle shall not be more than 15min;The maximum power limit is 50W.

**It is FORBIDDEN to overpower the antenna.**

Single antenna liver ablation is generally suitable for lesions **within 3cm**; If the lesion is larger than 3cm, double antenna ablation is recommended.





Ex Vivo Porcine Liver Ablation Data  
14G MRI Antenna on Continuous Mode with 4m cable



Liver

	50W(mm)	60W(mm)	80W(mm)	100W(mm)
3mins	<p>16mm 18mm</p>	<p>19mm 21mm</p>	<p>20mm 25mm</p>	<p>22mm 28mm</p>
5mins	<p>20mm 22mm</p>	<p>21mm 25mm</p>	<p>29mm 34mm</p>	<p>2.5mm 33mm</p>
8mins	<p>25mm 21mm</p>	<p>22mm 32mm</p>	<p>29mm 35mm</p>	<p>31mm 39mm</p>
10mins	<p>30mm 22mm</p>	<p>23mm 32mm</p>	<p>32mm 36mm</p>	<p>33mm 39mm</p>

**Note:**

The maximum power (100W) shown in the table is the limit power of the ablation antenna, which is not recommended to be used in clinical practice.





# Ex Vivo Porcine Liver Ablation Data 15G MRI Antenna on Continuous Mode with 4m cable



Liver

	30W(mm)	40W(mm)	50W(mm)
1min	<p>17mm 20mm</p>	<p>18mm 23mm</p>	<p>22mm 25mm</p>
3min	<p>18mm 21mm</p>	<p>20mm 25mm</p>	<p>24mm 28mm</p>
5min	<p>19mm 22mm</p>	<p>21mm 26mm</p>	<p>25mm 30mm</p>
8min	<p>22mm 25mm</p>	<p>25mm 30mm</p>	<p>30mm 37mm</p>
10min	<p>24mm 26mm</p>	<p>27mm 33mm</p>	<p>30mm 38mm</p>
12min	<p>25mm 33mm</p>	<p>30mm 35mm</p>	<p>31mm 40mm</p>
15min	<p>26mm 34mm</p>	<p>35mm 36mm</p>	<p>36mm 44mm</p>





Ex Vivo Porcine Liver Ablation Data  
15G MRI Antenna on Pulse Mode with 4m cable



Liver

	40W(mm)	60W(mm)	80W(mm)
1min	<p>7mm 12mm</p>	<p>10mm 20mm</p>	<p>13mm 21mm</p>
3min	<p>13mm 20mm</p>	<p>20mm 27mm</p>	<p>21mm 27mm</p>
5min	<p>20mm 20mm</p>	<p>25mm 28mm</p>	<p>26mm 32mm</p>
8min	<p>23mm 23mm</p>	<p>25mm 30mm</p>	<p>32mm 36mm</p>
10min	<p>24mm 25mm</p>	<p>30mm 32mm</p>	<p>35mm 38mm</p>
12min	<p>25mm 26mm</p>	<p>35mm 35mm</p>	<p>36mm 40mm</p>
15min	<p>27mm 33mm</p>	<p>36mm 38mm</p>	<p>38mm 45mm</p>






Ex Vivo Porcine Lung Ablation Data  
14G Ceramic Antenna on Continuous Mode (for lesion < 30mm)

	30W(mm)	40W(mm)	60W(mm)	80W(mm)
1min	<p>7mm 13mm</p>	<p>9mm 15mm</p>	<p>12mm 19mm</p>	<p>9mm 21mm</p>
3min	<p>15mm 17mm</p>	<p>12mm 25mm</p>	<p>17mm 24mm</p>	<p>18mm 35mm</p>
5min	<p>17mm 20mm</p>	<p>15mm 25mm</p>	<p>23mm 35mm</p>	<p>21mm 40mm</p>
8min	<p>18mm 22mm</p>	<p>21mm 33mm</p>	<p>23mm 37mm</p>	<p>23mm 42mm</p>
10min	<p>18mm 25mm</p>	<p>21mm 34mm</p>		
12min	<p>19mm 27mm</p>	<p>25mm 37mm</p>		





Lung



## Note:

Output power 50W~80W is routinely used in operation, and the continuous working time of a single ablation circle shall not be more than 15min;The maximum power limit is 80W.

**It is FORBIDDEN to overpower the antenna.**

Single antenna liver ablation is generally suitable for lesions **within 3cm**; If the lesion is larger than 3cm, double antenna ablation is recommended.





# Ex Vivo Porcine Lung Ablation Data 14G Ceramic Antenna on Pulse Mode (for lesion < 30mm)



Lung

	30W(mm)	40W(mm)	60W(mm)	80W(mm)
1min	<p>4mm 15mm</p>	<p>7mm 17mm</p>	<p>8mm 18mm</p>	<p>9mm 21mm</p>
3min	<p>8mm 20mm</p>	<p>12mm 25mm</p>	<p>16mm 30mm</p>	<p>17mm 32mm</p>
5min	<p>9mm 20mm</p>	<p>10mm 25mm</p>	<p>19mm 32mm</p>	<p>22mm 35mm</p>
8min	<p>10mm 21mm</p>	<p>15mm 30mm</p>	<p>20mm 34mm</p>	<p>23mm 37mm</p>
10min	<p>20mm 22mm</p>	<p>21mm 30mm</p>	<p>22mm 34mm</p>	<p>23mm 45mm</p>
12min	<p>21mm 26mm</p>	<p>22mm 33mm</p>	<p>23mm 35mm</p>	<p>28mm 47mm</p>





# Ex Vivo Porcine Lung Ablation Data 16G Ceramic Antenna on Continuous Mode



Lung

	20W(mm)	30W(mm)	40W(mm)	50W(mm)
1min	<p>4mm 11mm</p>	<p>5mm 18mm</p>	<p>7mm 20mm</p>	<p>8mm 21mm</p>
3min	<p>6mm 15mm</p>	<p>10mm 22mm</p>	<p>11mm 23mm</p>	<p>12mm 25mm</p>
5min	<p>7mm 18mm</p>	<p>15mm 24mm</p>	<p>18mm 30mm</p>	<p>19mm 31mm</p>
8min	<p>12mm 20mm</p>	<p>16mm 28mm</p>	<p>20mm 32mm</p>	<p>21mm 33mm</p>
10min	<p>13mm 21mm</p>	<p>17mm 30mm</p>	<p>23mm 34mm</p>	<p>26mm 36mm</p>

## Note:

Output power 30W~40W is routinely used in operation, and the continuous working time of a single ablation circle shall not be more than 15min;The maximum power limit is 50W.

**It is FORBIDDEN to overpower the antenna.**

Single antenna liver ablation is generally suitable for lesions **within 3cm**; If the lesion is larger than 3cm, double antenna ablation is recommended.





# Ex Vivo Porcine Lung Ablation Data 16G Ceramic Antenna on Pulse Mode



Lung

	20W(mm)	30W(mm)	40W(mm)	50W(mm)
1min	3mm  7mm	5mm  14mm	6mm  20mm	7mm  21mm
3min	4mm  7mm	6mm  18mm	7mm  21mm	8mm  22mm
5min	5mm  17mm	7mm  20mm	13mm  25mm	14mm  26mm
8min	6mm  20mm	15mm  21mm	18mm  26mm	19mm  28mm
10min	9mm  21mm	15mm  22mm	22mm  28mm	23mm  30mm

## Note:

Output power 30W~40W is routinely used in operation, and the continuous working time of a single ablation circle shall not be more than 15min;The maximum power limit is 50W.

It is **FORBIDDEN** to overpower the antenna.

Single antenna liver ablation is generally suitable for lesions **within 3cm**; If the lesion is larger than 3cm, double antenna ablation is recommended.





# Ex Vivo Porcine Liver Ablation Data 16G Metallic Antenna on Ablation Mode II



Thyroid

	30W(mm)	35W(mm)	40W(mm)
10s	<p>5mm 8mm</p>	<p>5mm 9mm</p>	<p>6mm 11mm</p>
15s	<p>5mm 9mm</p>	<p>6mm 10mm</p>	<p>7mm 13mm</p>
20s	<p>6mm 9mm</p>	<p>7mm 11mm</p>	<p>8mm 14mm</p>
25s	<p>7mm 9mm</p>	<p>8mm 11mm</p>	<p>9mm 15mm</p>

## Note:

Output power 30W~40W is routinely used in operation, and the continuous working time of a single ablation circle shall not be more than 25s;The maximum power limit is 40W. It is **FORBIDDEN to overpower the antenna.**

Single antenna thyroid ablation is generally suitable for lesions within 3cm; If the lesion is larger than 3cm, double antenna ablation is recommended.





# Ex Vivo Porcine Lung Ablation Data 17G Metallic Antenna on Ablation Mode II



Thyroid

	20W(mm)	25W(mm)	30W(mm)
10s	<p>4mm 7mm</p>	<p>4mm 8mm</p>	<p>5mm 9mm</p>
15s	<p>5mm 8mm</p>	<p>5mm 10mm</p>	<p>5mm 12mm</p>
20s	<p>6mm 11mm</p>	<p>6mm 13mm</p>	<p>7mm 15mm</p>
25s	<p>7mm 13mm</p>	<p>7mm 15mm</p>	<p>8mm 16mm</p>

## Note:

Output power 20W~25W is routinely used in operation, and the continuous working time of a single ablation circle shall not be more than 25s; The maximum power limit is 30W. It is **FORBIDDEN to overpower the antenna.**

Single antenna thyroid ablation is generally suitable for lesions within 3cm; If the lesion is larger than 3cm, double antenna ablation is recommended.





# Ex Vivo Porcine Liver Ablation Data 14G Hisphere Antenna on Continuous Mode

	1min	3min	5min	6min	7min	8min	10min	12min	15min
30W	 11mm (top), 11mm (bottom)	 14mm (top), 21mm (bottom)	 19mm (top), 22mm (bottom)	 20mm (top), 25mm (bottom)	 21mm (top), 27mm (bottom)	 22mm (top), 27mm (bottom)	 28mm (top), 32mm (bottom)	 30mm (top), 32mm (bottom)	 30mm (top), 33mm (bottom)
40W	 11mm (top), 19mm (bottom)	 16mm (top), 25mm (bottom)	 20mm (top), 28mm (bottom)	 24mm (top), 31mm (bottom)	 26mm (top), 31mm (bottom)	 30mm (top), 31mm (bottom)	 32mm (top), 33mm (bottom)	 36mm (top), 37mm (bottom)	 37mm (top), 37mm (bottom)
50W	 13mm (top), 20mm (bottom)	 21mm (top), 26mm (bottom)	 23mm (top), 28mm (bottom)	 26mm (top), 31mm (bottom)	 30mm (top), 31mm (bottom)	 32mm (top), 35mm (bottom)	 35mm (top), 36mm (bottom)	 38mm (top), 38mm (bottom)	 39mm (top), 40mm (bottom)
60W	 14mm (top), 24mm (bottom)	 21mm (top), 28mm (bottom)	 26mm (top), 30mm (bottom)	 32mm (top), 35mm (bottom)	 35mm (top), 36mm (bottom)	 35mm (top), 37mm (bottom)	 39mm (top), 40mm (bottom)	 39mm (top), 40mm (bottom)	 45mm (top), 47mm (bottom)
70W	 14mm (top), 24mm (bottom)	 21mm (top), 28mm (bottom)	 26mm (top), 32mm (bottom)	 32mm (top), 36mm (bottom)	 35mm (top), 37mm (bottom)	 35mm (top), 39mm (bottom)	 40mm (top), 40mm (bottom)	 41mm (top), 41mm (bottom)	 45mm (top), 47mm (bottom)
80W	 14mm (top), 26mm (bottom)	 26mm (top), 34mm (bottom)	 29mm (top), 34mm (bottom)	 32mm (top), 36mm (bottom)	 35mm (top), 38mm (bottom)	 35mm (top), 39mm (bottom)	 40mm (top), 40mm (bottom)	 41mm (top), 41mm (bottom)	 47mm (top), 49mm (bottom)
100W	 19mm (top), 26mm (bottom)	 27mm (top), 34mm (bottom)	 30mm (top), 36mm (bottom)	 32mm (top), 39mm (bottom)	 35mm (top), 40mm (bottom)	 36mm (top), 40mm (bottom)	 40mm (top), 44mm (bottom)	 41mm (top), 50mm (bottom)	 51mm (top), 54mm (bottom)

## Note:

Above data was obtained from ex-vivo porcine liver ablation experiment at the following cooling medium temperature conditions:  
Output power 30-60w-20°C room temperature saline water; Output power 70-100w-10°C cooled saline water.



**ECO MEDICAL TECHNOLOGY (NANJING) CO., LTD.**  
**(former: NANJING ECO MICROWAVE SYSTEM CO., LTD.)**

**Address:** 3rd & 4th Floors, J5 Building, NJUT Science & Technology Industrial Park, No.15 Wanshou Road, Jiangbei New Area, Nanjing, Jiangsu, China, 211800.

ECO Medical Instruments Industrial Park, No.19 Xinghui Road, Jiangbei New Area, Nanjing, Jiangsu, China, 210044.

**Tel:** +86-25-86262828

**Fax:** +86-25-86262777

**Email:** info@njeco.com.cn

**Website:** www.ecomicrowave.com



LinkedIn



ECO Website



Facebook



Youtube

Although we have tried our best to provide accurate and up-to-date information, no responsibility can be accepted for any mistakes contained in this catalogue. The products and descriptions, as shown in this catalogue may be changed without prior notice.