




ABLE[®] Haemodialysis Catheter kit



Focus on Quality, Respect for Life!

Haemodialysis Catheter



Soft Tip

Tapered tip facilitates ease of insertion without use of peel-away sheath. Reduces vessel trauma during insertion.



Side Holes

2.5cm far from end hole reduces recirculating and wall-suction. Minimizes risk of clot formation.



Rotatable Suture Wing







Facilitates skin inspection. Avoid blood vessel wall-suction. Minimizes risk of exit infections. No tube-kinking occurred.



Silicone Extension Tube

Increases patient comfort, will not crimp over time, able to visualize fluids.

Completeness

lumen	Pre-curved	Acute	Introducer Needle
 Single	Pre-curved	Acute	Introducer Needle
 Double I			
 Triple	Silicone extension tube	Chronic	Guide-wire

Specification

Product Code	Specification	Guidewire	Dilator	Introducer Needle
FH-1713	Single lumen 7F 10cm	0.035"×50cm	8F×10cm	18Ga
FH-1714	Single lumen 7F 13cm	0.035"×50cm	8F×10cm	18Ga
FH/FR-172913-5	Single lumen 7F 13.5cm	0.035"×50cm	8F×10cm	18Ga
FH-171915	Single lumen 7F 15cm	0.035"×50cm	8F×10cm	18Ga
FH-1715	Single lumen 7F 16cm	0.035"×50cm	8F×10cm	18Ga
FH-1716	Single lumen 7F 20cm	0.035"×60cm	8F×10cm	18Ga
FH-1813	Single lumen 8F 10cm	0.038"×50cm	9F×10cm	17Ga
FH-1814	Single lumen 8F 13cm	0.038"×50cm	9F×10cm	17Ga
FH/FR-182913-5	Single lumen 8F 13.5cm	0.038"×50cm	9F×10cm	17Ga
FH-181915	Single lumen 8F 15cm	0.038"×50cm	9F×10cm	17Ga
FH-1815	Single lumen 8F 16cm	0.038"×50cm	9F×10cm	17Ga
FH-1816	Single lumen 8F 20cm	0.038"×60cm	9F×10cm	17Ga
FH/FR-2611	Double lumen 6.5F 5cm	0.021"×50cm	7F×10cm	20Ga
FR-2611W	Double lumen 6.5F 5cm curved	0.021"×50cm	7F×10cm	20Ga
FH/FR-2612	Double lumen 6.5F 8cm	0.021"×50cm	7F×10cm	20Ga
FR-2612W	Double lumen 6.5F 8cm curved	0.021"×50cm	7F×10cm	20Ga
FH/FR-2613	Double lumen 6.5F 10cm	0.021"×50cm	7F×10cm	20Ga
FR-2613W	Double lumen 6.5F 10cm curved	0.021"×50cm	7F×10cm	20Ga
FH/FR-2614	Double lumen 6.5F 13cm	0.021"×50cm	7F×10cm	20Ga
FR-2614W	Double lumen 6.5F 13cm curved	0.021"×50cm	7F×10cm	20Ga
FH/FR-2615	Double lumen 6.5F 16cm	0.021"×50cm	7F×10cm	20Ga
FR-2615W	Double lumen 6.5F 16cm curved	0.021"×50cm	7F×10cm	20Ga
FH/FR-2911	Double lumen 8.5F 5cm	0.025"×50cm	9F×10cm	19Ga
FR-2911W	Double lumen 8.5F 5cm curved	0.025"×50cm	9F×10cm	19Ga
FH/FR-2912	Double lumen 8.5F 8cm	0.025"×50cm	9F×10cm	19Ga
FR-2912W	Double lumen 8.5F 8cm curved	0.025"×50cm	9F×10cm	19Ga
FH/FR-2913	Double lumen 8.5F 10cm	0.025"×50cm	9F×10cm	19Ga
FR-2913W	Double lumen 8.5F 10cm curved	0.025"×50cm	9F×10cm	19Ga
FH/FR-212911	Double lumen 8.5F 11cm	0.038"×60cm	9F×10cm	17Ga
FH/FR-212911W	Double lumen 8.5F 11cm curved	0.038"×60cm	9F×10cm	17Ga
FR-291912	Double lumen 8.5F 12cm	0.025"×50cm	9F×10cm	19Ga
FR-291912W	Double lumen 8.5F 12cm curved	0.025"×50cm	9F×10cm	19Ga
FH/FR-2914	Double lumen 8.5F 13cm	0.025"×50cm	9F×10cm	19Ga
FR-2914W	Double lumen 8.5F 13cm curved	0.025"×50cm	9F×10cm	19Ga
FR-291915	Double lumen 8.5F 15cm	0.025"×50cm	9F×10cm	19Ga
FR-291915W	Double lumen 8.5F 15cm curved	0.025"×50cm	9F×10cm	19Ga
FH/FR-2915	Double lumen 8.5F 16cm	0.025"×50cm	9F×10cm	19Ga
FR-2915W	Double lumen 8.5F 16cm curved	0.025"×50cm	9F×10cm	19Ga
FR-291917	Double lumen 8.5F 17cm	0.025"×50cm	9F×10cm	19Ga
FR-291917W	Double lumen 8.5F 17cm curved	0.025"×50cm	9F×10cm	19Ga
FR-2926	Double lumen 8.5F 20cm	0.025"×50cm	9F×10cm	19Ga
FR-2316	Double lumen 10F 20cm	0.038"×60cm	10F×15cm 9F×10cm	17Ga
FR-2313	Double lumen 10F 10cm	0.038"×50cm	10F×15cm 9F×10cm	17Ga
FR-2313W	Double lumen 10F 10cm curved	0.038"×50cm	10F×15cm 9F×10cm	17Ga
FR-2312	Double lumen 10F 8cm	0.038"×50cm	10F×15cm 9F×10cm	17Ga
FR-2312W	Double lumen 10F 8cm curved	0.038"×50cm	10F×15cm 9F×10cm	17Ga
FR-231912	Double lumen 10F 12cm	0.038"×50cm	10F×15cm 9F×10cm	17Ga
FR-231912W	Double lumen 10F 12cm curved	0.038"×50cm	10F×15cm 9F×10cm	17Ga

Specification

Product Code	Specification	Guidewire	Dilator	Introducer Needle
FR-231915	Double lumen 10F 15cm	0.038"×50cm	10F×15cm 9F×10cm	17Ga
FR-231915W	Double lumen 10F 15cm curved	0.038"×50cm	10F×15cm 9F×10cm	17Ga
FR-2314	Double lumen 10F 13cm	0.038"×50cm	10F×15cm 9F×10cm	17Ga
FR-2314W	Double lumen 10F 13cm curved	0.038"×50cm	10F×15cm 9F×10cm	17Ga
FR-231917	Double lumen 10F 17cm	0.038"×50cm	10F×15cm 9F×10cm	17Ga
FR-231917W	Double lumen 10F 17cm curved	0.038"×50cm	10F×15cm 9F×10cm	17Ga
FR-2315	Double lumen 10F 16cm	0.038"×50cm	10F×15cm 9F×10cm	17Ga
FR-2315W	Double lumen 10F 16cm curved	0.038"×50cm	10F×15cm 9F×10cm	17Ga
FR-2113	Double lumen 11.5F 10cm	0.038"×50cm	12F×15cm 10F×15cm	17Ga
FR-2113W	Double lumen 11.5F 10cm curved	0.038"×50cm	11.5F×15cm	17Ga
FR-212912	Double lumen 11.5F 12cm	0.038"×50cm	10F×15cm	17Ga
FR-212912W	Double lumen 11.5F 12cm curved	0.038"×50cm	12F×15cm 10F×15cm	17Ga
FR-211912G	Double lumen 11.5F 12cm pre-curved	0.038"×50cm	12F×15cm 10F×15cm	17Ga
FH/FR-2114	Double lumen 11.5F 13cm	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FR-2114W	Double lumen 11.5F 13cm curved	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FH/FR-212913-5	Double lumen 11.5F 13.5cm	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FH/FR-212913-5W	Double lumen 11.5F 13.5cm curved	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FH/FR-211915	Double lumen 11.5F 15cm	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FH/FR-211915W	Double lumen 11.5F 15cm curved	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FR-211915G	Double lumen 11.5F 15cm pre-curved	0.038"×50cm	12F×15cm 10F×15cm	17Ga
FH/FR-2115	Double lumen 11.5F 16cm	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FR-2115W	Double lumen 11.5F 16cm curved	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FR-2125W	Double lumen 11.5F 16cm curved	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FH/FR-2116	Double lumen 11.5F 20cm	0.038"×60cm	12F×15cm 10F×15cm	17Ga
FR-2116W	Double lumen 11.5F 20cm curved	0.038"×60cm	12F×15cm 10F×15cm	17Ga
FR-2115G	Double lumen 11.5F 20cm pre-curved	0.038"×60cm	12F×15cm 10F×15cm	17Ga
FR-211924	Double lumen 11.5F 24cm	0.038"×60cm	12F×15cm 10F×15cm	17Ga
FR-211924W	Double lumen 11.5F 24cm curved	0.038"×60cm	12F×15cm 10F×15cm	17Ga
FH/FR-2127	Double lumen 11.5F 25cm	0.038"×60cm	11.5F×15cm 10F×15cm	17Ga
FR-2117W	Double lumen 11.5F 25cm curved	0.038"×60cm	12F×15cm 10F×15cm	17Ga
FR-2213	Double lumen 12F 10cm	0.038"×50cm	12F×15cm 10F×15cm	17Ga
FR-2213W	Double lumen 12F 10cm curved	0.038"×50cm	12F×15cm 10F×15cm	17Ga
FR-221912G	Double lumen 12F 12cm pre-curved	0.038"×60cm	12F×15cm 10F×15cm	17Ga
FR-2216G	Double lumen 12F 20cm pre-curved	0.038"×60cm	12F×15cm 10F×15cm	17Ga
FH/FR-2214	Double lumen 12F 13cm	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FR-2214W	Double lumen 12F 13cm curved	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FH/FR-222913-5	Double lumen 12F 13.5cm	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FH/FR-222913-5W	Double lumen 12F 13.5cm curved	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FH/FR-221915	Double lumen 12F 15cm	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FH/FR-221915W	Double lumen 12F 15cm curved	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FR-221915G	Double lumen 12F 15cm pre-curved	0.038"×50cm	12F×15cm 10F×15cm	17Ga
FH/FR-2215	Double lumen 12F 16cm	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FR-2215W	Double lumen 12F 16cm curved	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FR-2225G	Double lumen 12F 16cm pre-curved	0.038"×50cm	12F×15cm 10F×15cm	17Ga
FG-2215	Double lumen 12F 16cm	0.038"×60cm	12F×15cm 10F×15cm	17Ga
FG-2215W	Double lumen 12F 16cm curved	0.038"×60cm	12F×15cm 10F×15cm	17Ga
FH/FR-2216	Double lumen 12F 20cm	0.038"×60cm	12F×15cm 10F×15cm	17Ga
FR-2216W	Double lumen 12F 20cm curved	0.038"×60cm	12F×15cm 10F×15cm	17Ga

Specification

Product Code	Specification	Guidewire	Dilator	Introducer Needle
FG-2216	Double lumen 12F 20cm	0.035"×60cm	12F×15cm 10F×15cm	17Ga
FG-2216W	Double lumen 12F 20cm curved	0.035"×60cm	12F×15cm 10F×15cm	17Ga
FR-221924	Double lumen 12F 24cm	0.035"×60cm	12F×15cm 10F×15cm	17Ga
FR-221924W	Double lumen 12F 24cm curved	0.035"×60cm	12F×15cm 10F×15cm	17Ga
FR-2217	Double lumen 12F 25cm	0.035"×60cm	12F×15cm 10F×15cm	17Ga
FR-2217W	Double lumen 12F 25cm curved	0.035"×60cm	12F×15cm 10F×15cm	17Ga
FR-2414	Double lumen 14F 13cm	0.038"×50cm	14F×15cm 12F×15cm	17Ga
FR-241915	Double lumen 14F 15cm	0.038"×60cm	14F×15cm 12F×15cm	17Ga
FR-241915W	Double lumen 14F 15cm curved	0.038"×60cm	14F×15cm 12F×15cm	17Ga
FR-2415	Double lumen 14F 16cm	0.038"×60cm	14F×15cm 12F×15cm	17Ga
FR-2415W	Double lumen 14F 16cm curved	0.038"×60cm	14F×15cm 12F×15cm	17Ga
FR-252919	Double lumen 14F 19cm	0.038"×60cm	11.5F×15cm 10F×15cm	17Ga
FR-2416	Double lumen 14F 20cm	0.038"×70cm	14F×15cm 12F×15cm	17Ga
FR-2416W	Double lumen 14F 20cm curved	0.038"×70cm	14F×15cm 12F×15cm	17Ga
FR-241924	Double lumen 14F 24cm	0.038"×70cm	14F×15cm 12F×15cm	17Ga
FR-241924W	Double lumen 14F 24cm curved	0.038"×70cm	14F×15cm 12F×15cm	17Ga
FR-2417	Double lumen 14F 25cm	0.038"×70cm	14F×15cm 12F×15cm	17Ga
FR-3214	Triple lumen 12F 13cm	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FR-3214W	Triple lumen 12F 13cm curved	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FR-3224	Triple lumen 12F 13cm	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FR-3224W	Triple lumen 12F 13cm curved	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FR-321915	Triple lumen 12F 15cm	0.038"×50cm	12F×15cm 10F×15cm	17Ga
FR-321915W	Triple lumen 12F 15cm curved	0.038"×50cm	12F×15cm 10F×15cm	17Ga
FR-3215	Triple lumen 12F 16cm	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FR-3215W	Triple lumen 12F 16cm curved	0.038"×50cm	11.5F×15cm 10F×15cm	17Ga
FR-3216	Triple lumen 12F 20cm	0.038"×60cm	12F×15cm 10F×15cm	17Ga
FR-3216W	Triple lumen 12F 20cm curved	0.038"×60cm	12F×15cm 10F×15cm	17Ga
FR-321924	Triple lumen 12F 24cm curved	0.038"×60cm	12F×15cm 10F×15cm	17Ga
FR-321924W	Triple lumen 12F 24cm	0.038"×60cm	12F×15cm 10F×15cm	17Ga
FR-3217W	Triple lumen 12F 25cm curved	0.038"×60cm	12F×15cm 10F×15cm	17Ga
FR-3217	Triple lumen 12F 25cm	0.038"×60cm	12F×15cm 10F×15cm	17Ga

Remarks:

1. The product codes in this chart are for HC kits with straight introducer needle configuration.
2. HC kits with Y-shaped introducer needle configuration are also available upon discussion. The product code of HC kit with Y-shaped introducer needle configuration is different from the product code of HC kit with straight introducer needle configuration. Before you order, please clarify which configuration you request, either straight or Y-shaped introducer needle.
3. Pediatric types of HC include double lumen 6.5Fr and 8.5Fr.
4. Adult types of HC include single lumen 7Fr, 8Fr; double lumen 10Fr, 11.5Fr, 12Fr, 14Fr and triple lumen 12Fr.
5. "FR" in the product code means with soft tip.
6. "FH" in the product code means without soft tip (i. e. with relatively harder tip).
7. Pre-curved catheter is available for adult type in double lumen 11.5Fr, 12Fr, 14Fr.



GUANGDONG BAIHE MEDICAL TECHNOLOGY CO.,LTD.

No. 89, Taoyuan East Road, Nanhai, Foshan, Guangdong Province, P.R.C.

Tel: +86-757-81207300 +86-757-81207302

Fax: +86-757-81207311

Email: able@baihemedical.com sales@baihemedical.com

<http://www.baihemedical.com>



CMC Medical Devices & Drugs S.L.
C/ Horacio Lengo No 18, CP 29006, Málaga, Spain.

