

EXTRACORPOREAL BLOOD CIRCUIT FOR BLOOD PURIFICATION EQUIPMENT



- Various models with high compatibility: Over ten models are compatible with most dialysis machines in the market.
- A wide range of accessories: Accessories such as infusion line, pre-fill connectors and waste bags are available for you to choose at will, and can meet the different demands.
- The products are certified by CFDA and EU authority.

Specifications:

Models	Pump tube&vein pot diameter parameters	Suitable equipment brand
JRHLS-001 Series	— —	Infusion tube for HDF/HF, tailored CRRT tube
JRHLL-010 Series	Pump tube:8*12mm Vein pot diameter:20mm	Nikkiso\Nipro\Toray,etc.
JRHLL-020 Series	Pump tube:8*12mm Vein pot diameter:22mm	Gambro\Fresenius\B. Braun, etc.
JRHLL-040 Series	Pump tube:8*12mm Vein pot diameter:30mm	Fresenius series

A.V. FISTULA NEEDLE SET

Ultra-thin dual radians sharp needle

Reduce tissue damage and pain

Rotation of the needle-wing and oval-shaped back hole

Adjustable to obtain the optimal blood flow and venous pressure,improve the quality of dialysis

Fine and even siliconized

Good biocompatibility, reduce puncture resistance

Red and blue points mark

Easy identification of needle point direction (red dot needle bevel down, blue dot needle bevel side up)



Technical Parameters

SPECIFICATIONS(LENGTH)	THE LENGTH OF CATHETER	CATEGORY	STERILIZATION METHOD	QUANTITY PER CARTON
15/16/17G (25mm)	300mm	Fixed wing needle	ETO	288/480pcs
15/16/17G (25mm)	300mm	Rotary wing needle	ETO	288/480pcs
15/16/17G (32mm)	300mm	Fixed wing needle	ETO	288/480pcs
15/16/17G (32mm)	300mm	Rotary wing needle	ETO	288/480pcs